

Lab Activity Handout **Answers**

Problem: Which nanoparticle will bleach (or “photo-sanitize”) water the fastest after UV light exposure: titanium dioxide, zinc oxide, or magnesium oxide?

Hypothesis: I predict that the titanium dioxide nanoparticle will bleach the fastest because the chemical properties of titanium dioxide make it the most effective.

Materials:



graduated cylinder



methyl orange



methylene blue



micro pipette




marker



stopwatch



plastic cups (2-oz)



distilled water

ZnO zinc oxide


MgO magnesium oxide

TiO₂ titanium dioxide

SAFETY
Equipment



goggles



lab coat/apron

Procedure:

1. Collect all materials ; designate responsibilities to each lab group member, if needed.
2. Obtain 8 plastic cups with lids and pour 20 ml of distilled water into each cup.
3. Pour 3 ml of methyl orange into the 4 cups with water.
4. Pour 3 ml of methylene blue into the other 4 cups with water.
5. Label the four methyl orange cups as ZnO _____, MgO _____, TiO₂ _____, & “CONTROL.”
6. Label the four methylene blue cups as ZnO _____, MgO _____, TiO₂ _____, & “CONTROL.”
7. Use a pipette to place 3 drops of each sample oxide as labeled on your cups & be sure to stir/mix your solutions well.
8. Take a picture of your methyl orange and methylene blue labeled cups with your phone (before light exposure) and be ready to take your cups outside for UV light exposure.
9. Using a stopwatch, record the time in seconds it takes for each sample to bleach (do not run longer than 10 min). Take another picture of the cups (after light exposure).

Data Table:

COMPLETE BLEACHING AFTER UV EXPOSURE in seconds				
Sample solutions	CONTROL	Titanium Dioxide	Magnesium Oxide	Zinc Oxide
Methylene blue	No bleaching			
Methyl orange	No bleaching			

Illustration: (students color in their results using their picture before & after light exposure.)

Samples before UV light exposure

Methylene Blue Samples



ZnO MgO TiO₂

Methyl Orange Samples



ZnO MgO TiO₂

Samples after UV light exposure

Methylene Blue Samples



ZnO MgO TiO₂

Methyl Orange Samples



ZnO MgO TiO₂

Conclusion: I accept my hypothesis because titanium dioxide was the nanoparticle that bleached the dye the fastest. This is due to its photocatalytic properties that we discussed in class.